

Automatic Autoclave TRA-603

## Automatic Autoclave TRA-603

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

Automatic Autoclave TRA-603 is a microprocessor-controlled tabletop electronic device with fully automatic, pre-programmed sterilization cycles and an LCD screen that displays continuous cycle information. It has an 18 L capacity, a maximum temperature of 134 °C, and a pressure of 210 MPa. It has an electro-polished stainless steel chamber and door. A simple user interface is very easy to use.

## Features:

- Tabletop autoclaves that are full auto and have pre-programmed sterilization cycles
- The chamber and entrance are made of stainless steel and have an electro-polish finish
- Control system with high precision for perfect sterilization results
- The built-in independent steam generator ensures that there is always steam accessible for faster cycles
- The user interface is very friendly and can operate very easily
- The LCD screen constantly provides data about the cycle
- The open-type water tank is very simple to fill and clean
- Microprocessor-controlled electronic autoclave that is fully automated
- It is a thermal autoclave with a non-vacuum working method
- Each operating status and program could be displayed as an icon on the LCD screen

## Specifications :

<b>Capacity</b>	18 L
<b>Dimension</b>	480×465×690 mm
<b>Heating Power</b>	2200 W
<b>Maximum Pressure</b>	210 MPa
<b>Maximum Temperature</b>	134 ?
<b>Maximum Timer</b>	0 to 55 min
<b>Number Of Trays</b>	3
<b>Power Supply</b>	220V 50Hz /110V60Hz
<b>Sterilization Room Size</b>	Ø247×350
<b>Weight</b>	47 kg



Northeast McClain Road Bentonville AR 72712, USA

Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)