



Balance TRB-603

www.axylab.com | info@axylab.com

Overview

Balance TRB-603 is a compact device with high stability; the sensor uses a unique locking device with ac/dc level, an LCD backlight display, and a weighing reaction speed range of 0 to 2000 g with a precision of 10 mg. It is intended to be used in conjunction with a printer or computer.

Features:

- It has ac/dc, level, LCD backlight display, fast weighing reaction time, and high stability
- The sensor employs a one-of-a-kind locking mechanism
- Counting function, percentage, automatic calibration, multiple calibration options, calibration weights
- With unit conversion, stability, memory function, animal weighing, full range gross weight reduction, and so on
- With an RS232 interface for connecting to a printer or computer

Specifications :

Accuracy	10 mg
Minimum Reading	0.01g
Pan Size	∅130 mm
Range	0 to 2000 g
Size	310×235×150 mm
Weight	1 Kg



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com