



## Centrifuge TRC-601

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

Centrifuge TRC-601 is a small device that can accommodate a variety of sample tubes. It has a maximum speed of 4000 rpm and an RCF of 1790× g. It has a rotor capacity of 20 ml×6. It is very easy to use, comes with a switch to continue operating and set parameters, and is well suited for qualitative analysis.

## Features:

- It is having with a Scope of supply that includes the rotor
- It can hold various sample tubes
- It is having with a Small footprint
- Very simple to operate
- Having a Step speed adjustment
- It Suits for qualitative analysis

## Specifications :

<b>Dimension</b>	Ø250×220 mm
<b>Maximum RCF</b>	1790×g
<b>Maximum RPM</b>	4000 RPM, step speed adjustment
<b>Rotor Capacity</b>	20 ml×6
<b>Timer</b>	0 to 30min
<b>Voltage</b>	110V 60Hz /220V 50Hz
<b>Weight</b>	3.5 Kg



Northeast McClain Road Bentonville AR 72712, USA  
Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)