



Centrifuge TRC-602

www.axylab.com | info@axylab.com

Overview

Centrifuge TRC-602 is a small device that can accommodate a variety of sample tubes. It has a maximum speed of 4000 rpm and an RCF of 1790× g. It has a rotor capacity of 20 ml×6. It is very easy to use, comes with a switch to continue operating and set parameters, and is well suited for qualitative analysis.

Features:

- It is having with a Scope of supply that includes the rotor
- It can hold various sample tubes
- It is having with a Small footprint
- Very simple to operate
- Having a Step speed adjustment
- It Suits for qualitative analysis

Specifications :

Dimension	Ø250×220 mm
Maximum RCF	1790×g
Maximum RPM	4000 RPM, step speed adjustment
Rotor Capacity	20 ml×6
Timer	0 to 30min
Voltage	110V 60Hz /220V 50Hz
Weight	3.5 Kg



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com