



Centrifuge TRC-603

www.axylab.com | info@axylab.com

Overview

Centrifuge TRC-603 with a powerful plastic sleeve and a steel lid, as well as a high brushless motor. It has a large LCD, a mechanical door lock, and a safety switch. It has a maximum rpm of 4000 and an RCF of 2325×g. It has an adjustable speed up and down and a standard rotor capacity of 20 ml ×12. It is appropriate for the qualitative analysis.

Features:

- It is a durable plastic container with a metal lid
- Having a brushless high-speed motor
- With an LED display, speed, and timer
- It has over speed and imbalance warnings
- It is built with a mechanical door lock and a safety switch are included
- It displays the error code and allows you to adapt the speed up and down
- It is appropriate for qualitative analysis

Specifications :

Dimension	280×310×265 mm
Maximum RCF	2325×g
Maximum RPM	4000 RPM, step-less speed adjustment, LED display speed, and timer
Rotor Capacity	20ml ×12 (standard)
Timer	0 to 99 min
Voltage	110 V 60 Hz /220 V 50 Hz
Weight	9 kg



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com