



Centrifuge TRC-604

www.axylab.com | info@axylab.com

Overview

Centrifuge TRC-604 is a brushless motor device with durable plastic housing and a metal lid. It has an LCD as well as an over-speed and imbalance alarm. It has a speed of 4000 RPM, an RCF of 2325 g, and a timer that runs from 0 to 99 minutes. It has a touch sensor control, and instantaneous stop when the lid is opened, and a qualitative analysis fits.

Features:

- It is a sturdy plastic case with a metal lid
- It has a brushless motor, an LCD for speed, RCF, and a timer, as well as an over-speed and imbalance alarm
- It comes with a touch panel control and an automatic stop when the lid is opened, as well as an error code display and speeds up or down adjustment
- suitable for qualitative analysis

Specifications :

Dimension	320×340×280 mm
Maximum RCF	2325×g
Maximum RPM	4000 RPM, step-less speed adjustment, LCD speed, RCF, and timer
Rotor Capacity	20ml×12 (standard)
Timer	0 to 99min
Voltage	110V 60Hz /220V50Hz
Weight	9kg



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com