



## **Fume Hood TRFH-603**

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

Fume Hood TRFH-603 is a ventilation cabinet that combines water, electricity, gas. The window is equipped with an automatic balance type. It is having with a 1.8M fume hood and having with a ventilation capacity is 450 to 1600 m<sup>3</sup>/h. A stainless steel water tank with liquid antiseptic on the bottom. Internal air discharges at a variable speed of 0.3m to 0.6m/s by the ventilation cabinet. A low-speed exhaust fan is included.

## Features:

- Water, electricity, gas, and air conditioning are all in the cabinet
- Inside, there is a multifunctional socket
- A stainless steel water tank with liquid antiseptic is installed on the bottom
- It is secure and dependable
- The ventilation capacity is 450 to 1600 m<sup>3</sup>/h, and the air circulation effectiveness is high
- It has a sliding glass door that can be opened and closed
- The flow plate makes bad gases at different heights discharged from different parts
- Using a water channel to solve a problem
- Comes with a 1.8 M fume cupboard
- It has a cabinet for chemical fumes
- It is an all-stainless steel fume cupboard

Specifications :

Dimension Of Working Area	1630×700×600 mm
Discharging Speed	within 0.3 m to 0.6 m/s adjustable
Fume Hood	1.8M
Inlet Pressure	>0.5 Pa
Overall Dimension W×D×H	1800×760×1950 mm
Power	800 W
Power Supply	AC, 220 V/50 Hz
Weight	200 Kg



Northeast McClain Road Bentonville AR 72712, USA  
Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)