



Digital Heating Water Bath TRHB-601

www.axylab.com | info@axylab.com

Overview

Digital Heating Water Bath TRHB-601 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 3 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

- It is equipped with a microprocessor control
- mechanism with a stainless steel interior
- Designed with a digital display inner temperature
- It is covered with an aluminum ring
- It is a thermostatic device

Specifications :

Capacity	3 L
Dimension	23x23x24 cm
Power	300 W
Temperature Fluctuation	±0.5 °
Temperature Homogeneity	±1 °
Temperature Range	(RT+5) -300 °
Voltage	110V 60Hz /220V 50 Hz
Weight	3.5kg
Working Size	16x16x12 cm



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com