



Digital Heating Water Bath TRHB-604

www.axylab.com | info@axylab.com

Overview

Digital Heating Water Bath TRHB-604 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °. It holds with a capacity of 17 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

- It is equipped with a microprocessor control
- mechanism with a stainless steel interior
- Designed with a digital display inner temperature
- It is covered with an aluminum ring
- It is a thermostatic device

Specifications :

Capacity	17 L
Dimension	53x37x24 cm
Power	1500 W
Temperature Fluctuation	±1.5 ?
Temperature Homogeneity	±1 ?
Temperature Range	(RT+5) -300 ?
Voltage	110V 60Hz /220V 50 Hz
Weight	7.5 kg
Working Size	46x30x12 cm



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com