



N Vertical Autoclave TRVA-601

www.axylab.com | info@axylab.com

Overview

N Vertical Autoclave TRVA-601 is made entirely of stainless steel, is simple to use, safe, and dependable, and comes with two stainless steel sterilizing baskets. Presented with a 35 L chamber volume, the working temperature was 134 °C, and the working pressure was 0.22 MPa. If no beeping reminder, the machine will shut down automatically. A digital display of the working status, and a touch key panel.

Features:

- Design made completely of impeccable high-combination steel
- Computerized show of working status, a dash of key
- After cleansing, the cool air and steam are naturally released
- After cleansing, the machine will naturally close down with a blare to remind you
- Outfitted with two treated steel disinfecting containers
- It is not tough to work, protected and reliable

Specifications :

Adjustment Of Temperature	105 to134 °C
Chamber Volume	35 L
G.W	56 Kg
Heat Average	?±1 °C
Maximum Working Pressure	0.23 MPa
N.W	42 Kg
Overall Dimension	480×460×850 mm
Power	AC220V.50Hz
Power Consumption	2.5Kw
Timer	0 to 99 min or 0 to 99 hour 59 min
Transport Dimension	570×550×970 mm
Working Pressure	0.22 MPa
Working Temperature	134 °C



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com