



## NC Linear Shaker TRLS-601

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

NC Linear Shaker TRLS-601 has dual LCD screens for parameter setting and viewing. It can withstand a maximum load weight and operates in a timed and constantly functioning mode. The dual LCD screens separately display time and speed, with a speed range of 40 to 200 rpm and a time range of 1 minute to 19 hours 59 minutes.

## Features:

- There is a maximum load-bearing weight of 7.5kg allowed
- The period is from 1min 19 hrs. 59 min
- The time and speed are displayed separately on dual LCD screens
- There are numerous inserts and fittings available
- Designed with a brushless DC motor requires no maintenance

Specifications :

Ambient Temperature	5 to 40°C
Amplitude	20mm
Humidity	0.8
Maximum Load	3kg
Motor Type	DC motor
Speed Display	LCD
Speed Range	40 to 200rpm
Time setting range	0 to 19h 59min
Voltage	100 to 240V, 50/60Hz



Northeast McClain Road Bentonville AR 72712, USA  
Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)