



NC Linear Shaker TRLS-602

www.axylab.com | info@axylab.com

Overview

NC Linear Shaker TRLS-602 has dual LCD screens for parameter setting and viewing. It can withstand a maximum load weight and operates in a timed and constantly functioning mode. The dual LCD screens separately display time and speed, with a speed range of 100 to 350 rpm and a time range of 1 minute to 19 hours 59 minutes. Present with the brushless DC motor is maintenance-free.

Features:

- There is a maximum load-bearing weight of 7.5kg allowed
- The period is from 1min 19 hrs. 59 min
- The time and speed are displayed separately on dual LCD screens
- There are numerous inserts and fittings available
- Designed with a brushless DC motor requires no maintenance

Specifications :

Ambient Temperature	5 to 40°C
Amplitude	10mm
Humidity	0.8
Maximum Load	7.5kg
Motor Type	DC motor
Speed Display	LCD
Speed Range	100 to 350rpm
Time setting range	0 to 19h 59min
Voltage	100 to 240V, 50/60Hz



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com