



## **Portable Sterilizer TRPS-602**

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

Portable Sterilizer TRPS-602 is a stainless steel device with an electric heating mode and a Dual scale numerical indication pressure gauge. It is having with a movable interior that is convenient for operation and with a capacity of 24 L and disinfecting Pressure of 0.14-0.165Mpa. It is having with the self-inflating type seal.

## Features:

- It is made of stainless steel
- It has an overpressure auto-discharging system. 0.145 to 0.165MPa
- The maximum working temperature is 126 °C
- Numeric indication with a dual magnitude Seal that self-inflates
- There are two modes: electric heating and electric/coal heating
- The movable interior makes operation easier

## Specifications :

<b>Capacity</b>	24 L
<b>Dimension</b>	39×39×65 cm
<b>Disinfecting Pressure</b>	0.14-0.165 MPa
<b>Heating Mode</b>	Electric-heating
<b>Maximum Working Temperature</b>	126 °
<b>Power</b>	2 KW
<b>Voltage</b>	110 V 60 Hz /220 V 50 Hz
<b>Weight</b>	19 kg



Northeast McClain Road Bentonville AR 72712, USA  
Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)