



Portable Sterilizer TRPS-603

www.axylab.com | info@axylab.com

Overview

Portable Sterilizer TRPS-603 is a stainless steel device with an electric/coal heating mode and a Dual scale numerical indication pressure gauge. It is having with a movable interior that is convenient for operation and with a capacity of 18 L and disinfecting Pressure of 0.14-0.165Mpa. It is having with the self-inflating type seal.

Features:

- It is made of stainless steel
- It has an overpressure auto-discharging system. 0.145 to 0.165MPa
- The maximum working temperature is 126 °C
- Numeric indication with a dual magnitude Seal that self-inflates
- There are two modes: electric heating and electric/coal heating
- The movable interior makes operation easier

Specifications :

Capacity	18 L
Dimension	39×39×44 cm
Disinfecting Pressure	0.14-0.165 MPa
Heating Mode	Electric/coal heating
Maximum Working Temperature	126 ?
Power	2 KW
Voltage	110 V 60 Hz /220 V 50 Hz
Weight	12.5 kg



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com