



Ultrasonic Cleaner TRUC-602

www.axylab.com | info@axylab.com

Overview

Ultrasonic Cleaner TRUC-602 has an ultrasonic vibration frequency of interaction composed of high-quality SUS304 stainless steel stampings. It has a 3 L capacity, a heating temperature range of 20 to 80°C, and a digital timer with a setting range of 0 to 99 minutes. It does not have a drainage ball valve.

Features:

- The intensity of ultrasonic vibration, thoroughly solve the ultrasonic cleaning
- Its work-performed tank is made of high-quality SUS304 stainless steel molding
- The bladder is made of high-quality steel material, and the cladding is made of aluminum
- Electronic timer with a period range of 0 to 99 minutes
- The drainage ball valves are installed in drainage systems 6L and higher models

Specifications :

| | |
|-----------------------------|---------------|
| 1 Unit Packing | 34×25×32 cm |
| 1 Unit Weight | 4 Kg |
| Dimensions | 27×25×32 Cm |
| Drainage | Without |
| Heating Power | 100W |
| Heating Temperature | 20 To 80° |
| Overall Packing | 79×38×69 Cm |
| Rinse Tank Size | 24×13.5×10 cm |
| Timer | 0 To 99min |
| Ultrasonic Frequency | 28/40khz |
| Ultrasonic Power | 120W |
| Volume | 3 L |
| Weight | 25Kg |



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com