



## Ultrasonic Cleaner TRUC-603

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

Ultrasonic Cleaner TRUC-603 has an ultrasonic vibration frequency of interaction composed of high-quality SUS304 stainless steel stampings. It has a 6 L capacity, a heating temperature range of 20 to 80°C, and a digital timer with a setting range of 0 to 99 minutes. It does not have a drainage ball valve.

## Features:

- The intensity of ultrasonic vibration, thoroughly solve the ultrasonic cleaning
- Its work-performed tank is made of high-quality SUS304 stainless steel molding
- The bladder is made of high-quality steel material, and the cladding is made of aluminum
- Electronic timer with a period range of 0 to 99 minutes
- The drainage ball valves are installed in drainage systems 6L and higher models

Specifications :

1 Unit Packing	43×28×42 cm
1 Unit Weight	6 Kg
Dimensions	33×18×31 Cm
Drainage	With
Heating Power	300W
Heating Temperature	20 To 80°
Overall Packing	86×46×82 Cm
Rinse Tank Size	30×15×15 cm
Timer	0 To 99min
Ultrasonic Frequency	28/40khz
Ultrasonic Power	180W
Volume	6 L
Weight	38 Kg



Northeast McClain Road Bentonville AR 72712, USA  
Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)