



## Water Distiller TRWD-601

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

Water Distiller TRWD-601 is corrosion and aging resistant and has an intelligent automatic control. It is equipped with a stainless steel bowl with an electric heater. It condenses during the steaming process for the distilled water by electric heating. It has an output water level of 5 liters per hour.

## Features:

- It is corrosion-resistant and age resistant
- Having a smart automatic control
- Electric heating causes steam condensation into distilled water
- It is outfitted with distillation generating equipment, which includes a condenser and other components
- It consists of a stain-free steel electric double distiller

Specifications :

Output	5 L/ hour
Package	300×300×630 mm
Power	AC220 V/ 4.5 KW
Weight	7 kg



Northeast McClain Road Bentonville AR 72712, USA  
Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)