



## **Water Distiller TRWD-603**

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

Water Distiller TRWD-603 is corrosion and aging resistant and has an intelligent automatic control. It is equipped with a stainless steel bowl with an electric heater. It condenses during the steaming process for the distilled water by electric heating. It has an output water level of 20 liters per hour.

## Features:

- It is corrosion-resistant and age resistant
- Having a smart automatic control
- Electric heating causes steam condensation into distilled water
- It is outfitted with distillation generating equipment, which includes a condenser and other components
- It consists of a stain-free steel electric double distiller

## Specifications :

<b>Output</b>	20L/hour
<b>Package</b>	350×495×790 mm
<b>Power</b>	AC380V/13.5 KW
<b>Weight</b>	10 kg



Northeast McClain Road Bentonville AR 72712, USA  
Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)