



Water Steam Distiller TRWD-604

www.axylab.com | info@axylab.com

Overview

Water Steam Distiller TRWD-604 is a stainless steel water tower system that is corrosion resistant and resistant to aging. Distillation of generating equipment, condenser, evaporator, electric heating pipe, and control unit. It has an automatic, intelligent control system. It is equipped with a control panel and has a button-type mode. Designed with a water level window, which can clearly show the potresidual water level.

Features:

- It is a water tower mechanism made from steel
- It is corrosion-resistant and age resistant
- It has an automatic and intelligent control system
- designed with a water level window that displays the residual water level
- The heating component is made up of immersion heating tubes that are long-lasting

Specifications :

Output	5 L/ hour
Package	540×270×490mm
Power	AC220 V/ 4.5 KW
Weight	9 Kg



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com