



## **Water Steam Distiller TRWD-605**

[www.axylab.com](http://www.axylab.com) | [info@axylab.com](mailto:info@axylab.com)

## Overview

Water Steam Distiller TRWD-605 is a corrosion-resistant as well as aging-resistant stainless steel mesh tower mechanism. Distillation of the generating equipment, including the condenser, evaporator, electric heating pipe, and control unit. It has an intelligent, automatic control system. It comes with a control panel and a button-type mode. Designed with a water level window that displays the residual water level.

## Features:

- It is a water tower mechanism made from steel
- It is corrosion-resistant and age resistant
- It has an automatic and intelligent control system
- designed with a water level window that displays the residual water level
- The heating component is made up of immersion heating tubes that are long-lasting

## Specifications :

<b>Output</b>	10 L/hour
<b>Package</b>	610×310×570mm
<b>Power</b>	AC380V/ 7.5 KW
<b>Weight</b>	11 Kg



Northeast McClain Road Bentonville AR 72712, USA  
Email: [info@axylab.com](mailto:info@axylab.com) | Website: [www.axylab.com](http://www.axylab.com)