



Water Steam Distiller TRWD-606

www.axylab.com | info@axylab.com

Overview

Water Steam Distiller TRWD-606 is a stainless steel water tower system that is resistant to corrosion and aging. Distillation of power plants, condensers, evaporators, electric heating pipes, and controls. It has an automatic and intelligent control system. It is equipped with a control panel and has a button mode. Equipped with a water level window that can display the remaining water level. The heating component is composed of an immersion heating tube and is durable.

Features:

- It is a water tower mechanism made from steel
- It is corrosion-resistant and age resistant
- It has an automatic and intelligent control system
- designed with a water level window that displays the residual water level
- The heating component is made up of immersion heating tubes that are long-lasting

Specifications :

Output	20L/hour
Package	740×400×680mm
Power	AC380V/13.5 KW
Weight	17 Kg



Northeast McClain Road Bentonville AR 72712, USA
Email: info@axylab.com | Website: www.axylab.com